

3008A Operation Table Hydraulic



TRIUP International Corp.

Floor#8, FengHuang Culture Plaza,
JiangDong Middle Road#215, Nanjing, P.R. China 210019
Tel: 0086 25 52334966 Fax: 0086 25 52334076
www.triup.com

Feature :

- Hydraulic control operation table with 4 section with kidney elevator.
- Mobile table with blocking system by a lever located in one side of the base.
- All table movement except legs section can be operated head end position by means two crank.
- Base cover and side rail made by stainless steel
- Table height adjustable hydraulically by means of foot pedal

Specification :

Dimension: 2100 ±50mm x 480±50mm (length x width)

Height Adjustment : 760-1010mm

Forward tilting : $\geq 55^\circ$

Back ward tilting : $\geq 35^\circ$

Lateral tilting (right /left) : $25^\circ /25^\circ$

Back rest adjustment : Up 85° / down 10°

Leg rest adjustment : Down 90°

Head section adjustment : Up to 90°

Kidney (waist) elevation : 120mm

Mattress thickness : 50mm

Stainless steel (ASI316/304) side bar/rails

Accessories (with clamp) :

1x Anesthetic screen

1x Pair solder support

1x Pair body support

1x Pair arm support

1x Pair legs support

1x Safety belt (patient body strap)

1x Kidney elevator set

1x Mattress

TRIUP International Corp.

Floor#8, FengHuang Culture Plaza,
JiangDong Middle Road#215, Nanjing, P.R. China 210019
Tel: 0086 25 52334966 Fax: 0086 25 52334076
www.triup.com